
Product Data Sheet

Product Name: 5,7-Dichlorokynurenic acid sodium salt

Cat. No.: GC14598

Chemical Properties

Cas. No. 1184986-70-6

Chemical Name sodium 5,7-dichloro-4-oxo-1,4-dihydroquinoline-2-carboxylate

SMILES C1C=CC(Cl)=C(C(C=C2C([O-])=O)=O)C(N2)=C1.[Na+]Formula $C_{10}H_4Cl_2NNaO_3$ M.Wt 280.04

Solubility <28mg/ml in 1eq. NaOH Storage Desiccate at RT

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA