
Product Data Sheet

Product Name: NPEC-caged-noradrenalin

Cat. No.: GC15704

Chemical Properties

Cas. No. 1257323-83-3

Chemical Name (E)-1-(2-nitrophenyl)ethyl hydrogen ((R)-2-(3,4-dihydroxyphenyl)-2-hydroxyethyl)carbonimidate

SMILES CC(O/C(O)=N/C[C@@](O)([H])C1=CC(O)=C(O)C=C1)([H])C2=CC=CC=C2N(=O)=OFormula C₁₇H₁₈N₂O₇ M.Wt 362.33

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA