

---

**Product Data Sheet**

---

Product Name: Neuromedin C (porcine)

Cat. No.: GC16376

**Chemical Properties**

Cas. No. 81608-30-2

Formula C<sub>50</sub>H<sub>73</sub>N<sub>17</sub>O<sub>11</sub>S

M.Wt 1120.29

Solubility Soluble to 1 mg/ml in sterile water

Storage Desiccate at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure 

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA