

---

**Product Data Sheet**

---

Product Name: ML 337  
Cat. No.: GC17095

**Chemical Properties**

Cas. No. 1443118-44-2

Chemical Name (R)-(2-fluoro-4-((4-methoxyphenyl)ethynyl)phenyl)(3-hydroxypiperidin-1-yl)methanone

SMILES COC1=CC=C(C=C1)C#CC2=CC(F)=C(C(N3CCC[C@](O)([H])C3)=O)C=C2

Formula  $C_{21}H_{20}FNO_3$  M.Wt 353.39

Solubility <35.34mg/ml in DMSO; <7.07mg/ml in ethanol Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA