
Product Data Sheet

Product Name: DAPK Substrate Peptide

Cat. No.: GC17650

Chemical Properties

Cas. No. 386769-53-5

Formula $C_{70}H_{115}N_{25}O_{17}$

M.Wt 1578.82

Solubility Soluble to 1 mg/ml in Water

Storage Desiccate at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Synthetic peptide substrate for death associated protein kinase (DAPK) ($K_M = 9 \mu M$).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA