
Product Data Sheet

Product Name: TB Dry, Sterile Solution

Cat. No.: GC17845

Chemical Properties

Cas. No.

Formula M.Wt

Solubility Storage store at room temperature

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

contains a sterile solution of TB Dry powder (AB01965) at 55 grams per liter. per liter: 12g tryptone, 24g yeast extract, 4g alternate carbon source 14.8g potassium phosphate

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA