

---

## Product Data Sheet

---

Product Name: YT Broth, 2X powder blend

Cat. No.: GC18115

### Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO      Storage      dry at room temperature

General tips      For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition      Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

for the growth of E.coli and other prokaryotes typical concentration: 31g/liter store dry at room temperature

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA