
Product Data Sheet

Product Name: 3-Pyrimidin-2-yl-Propionic Acid

Cat. No.: GC18330

Chemical Properties

Cas. No. 439108-20-0

Chemical Name 2-pyrimidinepropanoic acid

SMILES OC(=O)CCc1ncccn1Formula $C_7H_8N_2O_2$

M.Wt 152.2

Solubility DMF: 20 mg/ml, DMSO: 20 mg/ml, PBS (pH 7.2): .25 mg/ml

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure [3-Pyrimidin-2-yl-Propionic Acid](#)**Background**

3-Pyrimidin-2-yl-Propionic Acid is a synthetic intermediate useful for pharmaceutical synthesis.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA