
Product Data Sheet

Product Name: Isoproterenol sulfate dihydrate

Cat. No.: GC19698

Chemical Properties

Cas. No. 299-95-6

SMILES OC1=CC(C(CNC(C)C)O)=CC=C1O.OC2=CC(C(CNC(C)C)O)=CC=C2O.OS(O)(=O)=OFormula C₂₂H₄₀N₂O₁₂S

M.Wt 556.62

Solubility insoluble in EtOH; >59.9 mg/mL in Water; >74.7 mg/mL in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA