
Product Data Sheet

Product Name: Poly(acrylic acid)

Cat. No.: GC19723

Chemical Properties

Cas. No. 9003-01-4

Formula (C₃H₄O₂)_n

M.Wt

Solubility

Storage

Store at RT

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA