

---

**Product Data Sheet**

---

Product Name: Triton x-100

Cat. No.: GC19778

**Chemical Properties**

Cas. No. 9002-93-1

Formula  $(C_2H_4O)_N \cdot C_{14}H_{22}O$ 

M.Wt

Solubility DMSO : 100 mg/mL (ultrasonic and warming and heat to 60°C)

Storage Store at RT

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure 

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA