
Product Data Sheet

Product Name: PAMAM dendrimer, ethylenediamine core, generation 1.0 solution
Cat. No.: GC20150

Chemical Properties

Cas. No. 142986-44-5

Formula $\text{NH}_2(\text{CH}_2)_2\text{NH}_2:(\text{G}=1); \text{dendriPAMAM}(\text{NH}_2)_8$ M.Wt 1429.85

Solubility Storage 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA