

Product Data Sheet

Product Name: Histone
 Cat. No.: GC20183

Chemical Properties

Cas. No. 9064-47-5

Formula

M.Wt

Solubility H₂O: soluble 10 mg/mL, clear to hazy, colorless to light yellow

Storage Store at 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

□□□ □□□□□□□□