

Product Data Sheet

Product Name: Protein A Agarose (Fast Flow)

Cat. No.: GC20211

Chemical Properties

Cas. No.

Formula M.Wt

Solubility Storage Store at 4°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Protein A Agarose (Fast Flow) (Immunoprecipitation, IP) (Co-IP)

Protein A Agarose mouse IgG_{2a}, IgG_{2b}, IgA, rabbit IgG, human IgG₁, IgG₂, IgG₄

Protein A 6% agarose beads (Fast Flow) 2ml Protein A Agarose 2.5mg Protein A 2 ml Protein A Agarose 20mg human IgG (Linear flow rate) 50-300cm/h

Protein A Agarose TBS 2ml 0.5ml Agarose beads

Protein A Agarose 100

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA