
Product Data Sheet

Product Name: 2-Picolinic acid methyl ester

Cat. No.: GC25013

Chemical Properties

Cas. No. 2459-07-6

Formula C₇H₇NO₂ M.Wt 137.14

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

2-Picolinic acid methyl ester (Methyl picolinate, 2-Pyridinecarboxylic acid methyl ester, Methyl pyridine-2-carboxylate) is used as pharmaceutical intermediate.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA