
Product Data Sheet

Product Name: 4,4'-Bis(4-aminophenoxy)biphenyl

Cat. No.: GC25017

Chemical Properties

Cas. No. 13080-85-8

Formula C₂₄H₂₀N₂O₂

M.Wt 368.43

Solubility DMSO: 74 mg/mL (200.85 mM);;

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

4,4'-Bis(4-aminophenoxy)biphenyl is a monomer for polyimide production.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA