
Product Data Sheet

Product Name: Iso-H7 dihydrochloride

Cat. No.: GC25530

Chemical Properties

Cas. No. 140663-38-3

Formula C₁₄H₁₉Cl₂N₃O₂S

M.Wt 364.29

Solubility DMSO: 73 mg/mL (200.39 mM);Water: 73 mg/mL (200.39 mM);Ethanol: 73 mg/mL (200.39 mM)

Store
Storage at -
20°CGeneral
tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Iso-H7 dihydrochloride is an inhibitor of phosphokinase C.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA