
Product Data Sheet

Product Name: Paroxetine mesylate

Cat. No.: GC25708

Chemical Properties

Cas. No. 217797-14-3

Formula C₁₉H₂₀FNO₃.CH₄O₃S

M.Wt 425.47

Solubility DMSO: 85 mg/mL (199.78 mM); Water: 10 mg/mL (23.50 mM);

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Paroxetine Mesylate (BRL-29060A mesylate, FG-7051 mesylate) is the mesylate salt form of paroxetine, a phenylpiperidine derivative and a selective serotonin reuptake inhibitor (SSRI) with antidepressant and anxiolytic properties.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA