
Product Data Sheet

Product Name: Piperazine phosphate

Cat. No.: GC25754

Chemical Properties

Cas. No. 85547-56-4

Formula C₂₉H₃₂Cl₂N₆.H₃O₄P

M.Wt 633.51

Solubility DMSO: Insoluble; Water: 2 mg/mL (3.16 mM); Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Piperazine phosphate is an orally active bisquinolone antimalarial drug.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA