

Product Data Sheet

Product Name: PKH 67 Cell Membrane Labeling Kit  
Cat. No.: GC26306

Components

Solution A PKH67 stock solution 0.1ml  
Solution B PKH67 dilution 10ml  
Storage -20°C, protect from light

Protocol

1. 细胞膜标记  
① 将Solution A (PKH67 stock solution) 加入含有细胞的培养液中，使终浓度为10<sup>7</sup>/ml。  
② 将培养液置于37°C培养5min，然后置于4°C培养15min。  
③ 离心400×g，5min。  
④ 将细胞重悬于100μl PBS中。  
⑤ 将细胞置于20~25°C培养15~30min。  
⑥ 离心400×g，10min，用PBS洗涤2-3次。  
⑦ 检测Ex/Em=490/502nm

Caution: Product has not been fully validated for medical applications. For research use only.  
Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com  
Address: 10292 Central Ave. #205, Montclair, CA, USA

