

---

**Product Data Sheet**

---

Product Name: CY3-YNE  
Cat. No.: GC30161

**Chemical Properties**

Cas. No. 1010386-62-5

SMILES CC(/C(N1CCCCC(NCC#C)=O)=C\C=C\C2=[N+](CC)C(C=CC(S(=O)(O)=O)=C3)=C3C2(C)C)(C)C4=C1C=CC(S(=O)(O)=O)=C4

Formula C<sub>34</sub>H<sub>42</sub>N<sub>3</sub>O<sub>7</sub>S<sub>2</sub> M.Wt 668.84

Solubility DMSO : ≥ 30 mg/mL (44.85 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

CY3-YNE is a dye for the labeling of soluble proteins, peptides, and oligonucleotides/DNA.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA