

---

**Product Data Sheet**

---

Product Name: D-Pipecolinic acid

Cat. No.: GC30354

**Chemical Properties**

Cas. No. 1723-00-8

SMILES O=C([C@@H]1NCCCC1)OFormula  $C_6H_{11}NO_2$ 

M.Wt 129.16

Solubility  $H_2O$  : 66.67 mg/mL (516.18 mM; Need ultrasonic); DMSO :  
< 1 mg/mL (insoluble or slightly soluble)Store  
Storage at -  
20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

D-Pipecolinic acid is a normal human metabolite found in human biofluids.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA