
Product Data Sheet

Product Name: D-Lys(Z)-Pro-Arg-pNA (Chromozym Pca)

Cat. No.: GC30441

Chemical Properties

Cas. No. 108963-69-5

SMILES {d-Lys(Z)}-Pro-Arg-pNA

Formula C₃₁H₄₃N₉O₇ M.Wt 653.73

Solubility Soluble in Water Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

D-Lys(Z)-Pro-Arg-pNA is a luminescent substrate of activated protein C (APC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA