
Product Data Sheet

Product Name: Sotirimod (R850)

Cat. No.: GC30463

Chemical Properties

Cas. No. 227318-75-4

SMILES NC1=NC2=CC=CN=C2C3=C1N=C(C)N3CC(C)CFormula $C_{14}H_{17}N_5$

M.Wt 255.32

Solubility DMSO : 25 mg/mL (97.92 mM; ultrasonic and warming and heat to 80°C)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Sotirimod is an immunostimulant, and can potentially treat for actinic keratosis.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA