
Product Data Sheet

Product Name: p38 MAPK-IN-2

Cat. No.: GC30609

Chemical Properties

Cas. No. 635725-16-5

SMILES O=C(N1CCC(C2=NNC(C3=CC=C(Cl)C=C3F)=C2C4=NC=NC=C4)CC1)COFormula C20H19ClFN5O2 M.Wt 415.85

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

p38 MAPK-IN-2 is an inhibitor of p38 kinase.

[1]. CRYSTALLINE FORM OF 2-(4-(4-CHLORO-2-FLUOROPHENYL)-4-PYRIMIDIN-4-YL-1H-PYRAZOL-5-YL)PIPERIDIN-1-YL-2-OXOETHANOL. WO/2005/061486 A1

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
