

---

**Product Data Sheet**

---

Product Name: Hosenkoside M ((+)-Hosenkoside M)

Cat. No.: GC30643

**Chemical Properties**

Cas. No.

Formula M.Wt

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Hosenkoside M is a baccharane glycoside isolated from the seeds of *impatiens balsamina*.

[1]. Yu X, et al. An LC-MS/MS method for simultaneous determination of hosenkoside A and hosenkoside K from *Semen Impatiensis* in rat plasma and its application to a pharmacokinetic study. *Biomed Chromatogr.* 2017 Mar;31(3).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA