
Product Data Sheet

Product Name: (5R)-BW-4030W92

Cat. No.: GC30701

Chemical Properties

Cas. No. 189013-61-4

SMILES NC1=NC(CF)=[C@]([C@]2=CC=CC(Cl)=C2Cl)C(N)=N1.[R]Formula $C_{11}H_9Cl_2FN_4$ M.Wt 287.12

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

(5R)-BW-4030W92 is the R enantiomer of BW-4030W92. BW-4030W92 is a non-selective, voltage-, and use-dependent sodium channel antagonist.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA