

---

**Product Data Sheet**

---

Product Name: Spirendolol (Li 32-468)

Cat. No.: GC31005

**Chemical Properties**

Cas. No. 81840-58-6

SMILES O=C1C2=C(C(OCC(O)CNC(C)(C)C)=CC=C2)CC31CCCCC3Formula  $C_{21}H_{31}NO_3$  M.Wt 345.48

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**Spirendolol is a  $\beta$  adrenergic receptor antagonist.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA