
Product Data Sheet

Product Name: Paliroden (SR 57667)

Cat. No.: GC31034

Chemical Properties

Cas. No. 188396-77-2

SMILES FC(C1=CC(C2=CCN(CCC3=CC=C(C4=CC=CC=C4)C=C3)CC2)=CC=C1)(F)FFormula C₂₆H₂₄F₃N M.Wt 407.47

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Paliroden is an orally bioactive neurotrophic, non-peptidic compound that activates synthesis of endogenous neurotrophines, used for treatment of Alzheimer's Disease and Parkinson's.

[1]. Labie C, et al. Effects of paliroden (SR57667B) and xaliproden on adult brain neurogenesis. *Curr Alzheimer Res.* 2006 Feb;3(1):35-6.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA