
Product Data Sheet

Product Name: LH secretion antagonist 1

Cat. No.: GC31226

Chemical Properties

Cas. No. 88531-67-3

SMILES OC12CCN(CCC#C)C(C2(C)C)CC3=CC=C(O)C=C13.[H]ClFormula $C_{18}H_{24}ClNO_2$ M.Wt 321.84

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

LH secretion antagonist 1 is an antagonist of luteinising hormone secretion, and may be used as an analgesic.

[1]. Method of inhibiting luteinizing hormone secretion with 6,7-benzomorphan derivatives. US4406904A

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
