

---

**Product Data Sheet**

---

Product Name: 5-HT3 antagonist 1

Cat. No.: GC31263

**Chemical Properties**

Cas. No. 129294-09-3

SMILES O=C(C1=NNC2=C1C=CC=C2)NC3CN(CC4=CC=CC(C)=C4)CCN(C)C3Formula C<sub>22</sub>H<sub>27</sub>N<sub>5</sub>O M.Wt 377.48

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

5-HT3 antagonist 1 is a potent and selective antagonist of serotonin 3 (5-HT3) receptor.

[1]. INDAZOLE-3-CARBOXYLIC ACID DERIVATIVES. US5017573

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA