

---

**Product Data Sheet**

---

Product Name: Piperylone (PR66)

Cat. No.: GC31297

**Chemical Properties**

Cas. No. 2531-04-6

SMILES O=C1N(C2CCN(C)CC2)NC(C3=CC=CC=C3)=C1CCFormula  $C_{17}H_{23}N_3O$  M.Wt 285.38

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Piperylone can be used as an antispasmodic.

[1]. Martínková J, et al. A study of the inhibition of adrenaline-induced vasoconstriction in the isolated perfused liver of rabbit. *Hepatology*. 1990 Nov;12(5):1157-65.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---