
Product Data Sheet

Product Name: GPR40 Agonist 2

Cat. No.: GC31634

Chemical Properties

Cas. No. 1147729-48-3

SMILES OC(C[C@H](C#CC)C(C=C1)=CC=C1OCC2=CCCC3(CCCCC3)C2)=OFormula C₂₄H₃₀O₃ M.Wt 366.49

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

GPR40 Agonist 2 is a GPR40 agonist that can be used in the research of diabetes, extracted from patent WO2009054479A1.

GPR40 Agonist 2 is a GPR40 agonist that can be used in the research of diabetes[1].

[1]. SHIMADA, Takashi. SPIRO-RING COMPOUND AND USE THEREOF FOR MEDICAL PURPOSES. WO 2009054479 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
