

---

**Product Data Sheet**

---

Product Name: Anti-inflammatory agent 2

Cat. No.: GC32009

**Chemical Properties**

Cas. No. 133012-00-7

SMILES O=C(NS(=O)(C)=O)C(C1=CC=C(OCC2=NC3=CC=CC=C3C=C2)C=C1)C4CCCCC4Formula  $C_{26}H_{30}N_2O_4S$  M.Wt 466.59

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Anti-inflammatory agent 2 can be used to treat and prevent inflammatory diseases extracted from patent WO 2001035936 A2, example 1.

[1]. William Taylor, et al. Use of substituted (quinolin-2-yl-methoxy)phenyl-acyl-sulphonamides and -cyanamides for the treatment of diseases. WO 2001035936 A2.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA