
Product Data Sheet

Product Name: Monobehenin

Cat. No.: GC32364

Chemical Properties

Cas. No. 30233-64-8

SMILES CCCCCCCCCCCCCCCCCCCC(OCC(O)CO)=OFormula $C_{25}H_{50}O_4$

M.Wt 414.66

Solubility DMSO : < 1 mg/mL (insoluble or slightly soluble); Water : < 0.1 mg/mL (insoluble)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Monodocosanoin is a monoacylglycerol that contains docosanoic acid at one position.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA