

---

## Product Data Sheet

---

Product Name: Pancreatin

Cat. No.: GC32412

### Chemical Properties

Cas. No. 8049-47-6

SMILES [Pancreatin]

Formula

M.Wt

Solubility DMSO : 5 mg/mL ;Water : < 0.1 mg/mL (insoluble) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Pancreatin (Pancrelipase, Pancreatic enzymes), are commercial mixtures of amylase, lipase, and protease. They are used to treat malabsorption syndrome due to pancreatic problems.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA