
Product Data Sheet

Product Name: Hepronicate (Megrin)

Cat. No.: GC32527

Chemical Properties

Cas. No. 7237-81-2

SMILES O=C(C1=CN=CC=C1)OCC(COC(C2=CN=CC=C2)=O)(CCCCCC)COC(C3=CN=CC=C3)=OFormula C₂₈H₃₁N₃O₆ M.Wt 505.56

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Hepronicate is a peripheral vasodilator with blood lipid lowering action.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA