
Product Data Sheet

Product Name: PRX933 hydrochloride (GW876167 hydrochloride)

Cat. No.: GC32595

Chemical Properties

Cas. No. 639029-42-8

SMILES C[C@H]1N(C2=NC=CN=C2OCCOC3=CC=CN=C3)CCNC1.ClFormula $C_{16}H_{22}ClN_5O_2$ M.Wt 351.83

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

PRX933 hydrochloride is a 5-HT_{2c} receptor agonist extracted from patent WO 2014140631 A1.

PRX933 hydrochloride is a 5-HT_{2c} receptor agonist which is useful for the acute treatment of hypertension.

[1]. Peter Richardson. Prx933 (prx-00933), a 5ht_{2c} agonist, for use in the treatment of hypertension. WO 2014140631 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
