

---

**Product Data Sheet**

---

Product Name: Fmoc-Val-Cit-PAB-PNP

Cat. No.: GC32800

**Chemical Properties**

Cas. No. 863971-53-3

SMILES Fmoc-Val-Cit-PAB-PNP

Formula C<sub>40</sub>H<sub>42</sub>N<sub>6</sub>O<sub>10</sub>

M.Wt 766.8

Solubility DMSO : ≥ 40 mg/mL (52.16 mM); Water : < 0.1 mg/mL  
(insoluble)Storage Store at -  
20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Val-Cit-PAB-PNP is a peptide prodrug linker, is a linker for antibody-drug-conjugation (ADC).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA