
Product Data Sheet

Product Name: E3 ligase Ligand-Linker Conjugates 6

Cat. No.: GC32839

Chemical Properties

Cas. No. 2097973-72-1

SMILES O=C([C@H]1N(C([C@@H](NC(COCCOCCN)=O)C(C)(C)C)=O)C[C@H](O)C1)NCC2=CC=C(C3=C(C)N=CS3)C=C2.[H]Cl

Formula C₂₈H₄₂ClN₅O₆S M.Wt 612.18

Solubility DMSO : 100 mg/mL (163.35 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

E3 ligase Ligand-Linker Conjugates 6 is a synthesized compound that incorporates an E3 ligase ligand and a linker used in PROTAC technology.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA