
Product Data Sheet

Product Name: E3 Ligase Ligand-Linker Conjugates 22

Cat. No.: GC32903

Chemical Properties

Cas. No. 2225940-52-1

SMILES O=C1N(C(CC2)C(NC2=O)=O)C(C3=C1C=CC=C3NCCOCCOCCOCCOCCN)=OFormula $C_{23}H_{32}N_4O_8$ M.Wt 492.52

Solubility DMSO : 50 mg/mL (101.52 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

E3 Ligase Ligand-Linker Conjugates 22 incorporates an E3 ligase ligand and a linker, can be used for the treatment of EZH2-mediated cancer.

[1]. ARYL SULFONOHYDRAZIDES. WO 2018081530 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA