
Product Data Sheet

Product Name: D-Lys(Z)-Pro-Arg-pNA diacetate (Spectrozyme PCa)

Cat. No.: GC33560

Chemical Properties

Cas. No. 108963-70-8

SMILES {d-Lys(Z)}-Pro-Arg-pNA

Formula C₃₅H₅₁N₉O₁₁

M.Wt 773.83

Solubility Water : ≥ 8.88 mg/mL (11.48 mM)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

D-Lys(Z)-Pro-Arg-pNA (diacetate) diacetate is a chromogenic substrate.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA