
Product Data Sheet

Product Name: Tyrosine Protein Kinase JAK 2 Phospho-Tyr8, 9

Cat. No.: GC33586

Chemical Properties

Cas. No. 247171-44-4

SMILES Val-Leu-Pro-Gln-Asp-Lys-Glu-{pTyr}{pTyr}-Lys-Val-Lys-Glu-Pro-Gly-Glu

Formula C₈₈H₁₃₈N₂₀O₃₄P₂ M.Wt 2082.1

Solubility Soluble in Water Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Tyrosine Protein Kinase JAK 2 (Phospho-Tyr8, 9) is a peptide corresponding to amino acids 475 to 491 of mouse JAK2.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA