
Product Data Sheet

Product Name: Diacylglycerol acyltransferase inhibitor-1

Cat. No.: GC33808

Chemical Properties

Cas. No. 1610800-25-3

SMILES O=C1N(C2=CC3=C(N(CCC)C=C3)C(C4CC4)=C2)CCOC5=NC=NC(N)=C15Formula $C_{21}H_{23}N_5O_2$ M.Wt 377.44

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Diacylglycerol acyltransferase inhibitor-1 is a diacylglycerol acyltransferase (DGAT1) inhibitor.

[1]. Diacylglycerol acyltransferase inhibitor-1 is a diacylglycerol acyltransferase (DGAT1) inhibitor.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA