

---

**Product Data Sheet**

---

Product Name: FR183998 free base

Cat. No.: GC33895

**Chemical Properties**

Cas. No. 239440-20-1

SMILES O=C(C1=CC(C2=C(CI)SC(CI)=C2)=CC(C(NCCN(C)C)=O)=C1)NC(N)=NFormula  $C_{17}H_{19}Cl_2N_5O_2S$  M.Wt 428.34

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Protocol****Cell experiment:**

Intracellular pH of rats lymphocytes is measured at 37°C with a spectrofluorometer by using the ratio of the emission (530 nm) obtained at 490- and 440-nm excitation wavelengths. After addition of 10 µL of the BCECF- and acid-loaded cells to 460 µL of Na-free buffer with 5 µL of DMSO or FR183998, 25 µL of 2.0 M NaCl (final 100 mM) is applied to start the reaction. The pHi change of cells is measured for >2 min. The initial increase in pHi in response to the added NaCl is taken as an estimate of Na<sup>+</sup>/H<sup>+</sup>-exchange activity, and the inhibitory effect of the drug is evaluated as IC<sub>50</sub> value. The BCECF fluorescence signals are calibrated by titration with stepwise addition of small volumes of 1 M MES after permeabilization of the cells with 0.5% Triton X-100[1].

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---

## Product Data Sheet

---

### Animal experiment:

Male Sprague-Dawley rats age 9-10 weeks are anesthetized with sodium pentobarbital (50 mg/kg, i.p.), intubated, and ventilated with room air at a stroke volume of 4.5 mL with a frequency of 50 strokes/min. Saline or FR183998 is injected into the left femoral vein (bolus i.v.) 5 min before occlusion. Thoracotomy for coronary artery ligation is performed at the left side of the thorax above the heart, and the origin of the left anterior descending coronary artery is occluded by applying negative tension by aspiration through polyethylene tubing connected to the vacuum system, which can achieve regional ischemia. Successful occlusion is confirmed by ischemia-induced alteration of electrocardiogram (ECG). After 5-min ischemia, reperfusion is performed, and the incidence of ventricular tachycardia (VT) and ventricular fibrillation (VF) of >5 s is observed intermittently over a 5-min period of reperfusion[1].

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---

## Product Data Sheet

---

### References:

[1]. Ohara F, et al. Preischemic and postischemic treatment with a new Na<sup>+</sup>/H<sup>+</sup>-exchange inhibitor, FR183998, shows cardioprotective effects in rats with cardiac ischemia and reperfusion. J Cardiovasc Pharmacol. 1999 Dec;34(6):848-56.

[2]. Ishizaki M, et al. Protective effect of FR183998, a Na<sup>+</sup>/H<sup>+</sup> exchanger inhibitor, and its inhibition of iNOS induction in hepatic ischemia-reperfusion injury in rats. Shock. 2008 Sep;30(3):311-7.

### Background

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---

## Product Data Sheet

---

### Background

FR183998 free base is a potent Na<sup>+</sup>/H<sup>+</sup>-exchange inhibitor, with IC<sub>50</sub>s of 0.3 nM, 3.1 nM and 6.5 nM by measurement of pHi change in rat lymphocytes, rat and human platelets, respectively.

FR183998 free base is a Na<sup>+</sup>/H<sup>+</sup>-exchange inhibitor, with IC<sub>50</sub>s of 0.3 nM, 6.5 nM and 3.1 nM by measurement of pHi change in rat lymphocytes, rat and human platelets, respectively[1].

FR183998 (0.1 and 1.0 mg/kg, i.v.) shows no effect hemodynamic parameters, and does not affect mean blood pressure and heart rate in conscious rats. Pretreatment of 0.01, 0.032, 0.10 mg/kg FR183998 or posttreatment of 0.032 and 0.10 mg/kg FR183998 via intravenous administration, dose-dependently reduces reperfusion-induced ventricular fibrillation (VF) and mortality in reperfusion-induced arrhythmias in anesthetized rats, with ED<sub>50</sub>s against VF of 0.015 mg/kg and 0.070 mg/kg, respectively. FR183998 also reduces myocardial infarct sizes, and suppresses the arrhythmias in anesthetized rats[1]. FR183998 (1 mg/kg, i.v.) reduces the increase in serum levels of alanine transaminase, aspartate transaminase, and lactate dehydrogenase induced by hepatic I/R, and prevents the incidences of hepatic necrosis, apoptosis, and neutrophil infiltration. FR183998 blocks the I/R-induced activation of the NF-κB, reduces induction of iNOS and inhibits the production of nitric oxide. FR183998 also decreases the expression of the iNOS gene antisense transcript in the liver of hepatic I/R rats[2].

[1]. Ohara F, et al. Preischemic and postischemic treatment with a new Na<sup>+</sup>/H<sup>+</sup>-exchange inhibitor, FR183998, shows cardioprotective effects in rats with cardiac ischemia and reperfusion. *J Cardiovasc Pharmacol.* 1999 Dec;34(6):848-56. [2]. Ishizaki M, et al. Protective effect of FR183998, a Na<sup>+</sup>/H<sup>+</sup> exchanger inhibitor, and its inhibition of iNOS induction in hepatic ischemia-reperfusion injury in rats. *Shock.* 2008 Sep;30(3):311-7.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA