
Product Data Sheet

Product Name: OT antagonist 3 analog

Cat. No.: GC34363

Chemical Properties

Cas. No. 2444044-37-3

SMILES CC1=NN=C(C2=CN=C(OC3=CC=CC(C)=C3C)C=N2)N1C(C=C4)=CNC4=OFormula $C_{20}H_{18}N_6O_2$ M.Wt 374.4

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

OT antagonist 3 analog is an analog of OT antagonist 3.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA