
Product Data Sheet

Product Name: E3 ligase Ligand-Linker Conjugates 15

Cat. No.: GC34365

Chemical Properties

Cas. No.

SMILES O=C(N(C(CCC1=O)C(N1)=O)C2=O)C(C2=CC=C3)=C3OCC(NCCCC(O)=O)=O

Formula C19H19N3O8 M.Wt 417.37

Solubility DMSO : ≥ 55.5 mg/mL (132.98 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

E3 ligase Ligand-Linker Conjugates 15 is a synthesized compound that incorporates an E3 ligase ligand and a linker used in PROTAC technology.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA