
Product Data Sheet

Product Name: Super-TDU 1-31

Cat. No.: GC34818

Chemical Properties

Cas. No.

Formula C₁₄₁H₂₁₈N₄₀O₄₈

M.Wt

3241.48

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Super-TDU (1-31) is a peptide of Super-TDU, which is an inhibitor of YAP-TEADs, shows potent anti-tumor activity[1].

[1]. Jiao S, et al. A peptide mimicking VGLL4 function acts as a YAP antagonist therapy against gastric cancer. Cancer Cell. 2014 Feb 10;25(2):166-80.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA