
Product Data Sheet

Product Name: T3-ATA S-isomer

Cat. No.: GC34822

Chemical Properties

Cas. No. 2438721-48-1

SMILES CC(SCC(N[C@H](C(O)=O)CC(C=C1I)=CC(I)=C1OC2=CC=C(O)C(I)=C2)=O)=O

Formula $C_{19}H_{16}I_3NO_6S$ M.Wt 767.11

Solubility Soluble in DMSO Storage 4°C, protect from light, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

T3-ATA S-isomer is the S-isomer of T3-ATA, which is the active form of the thyroid hormone.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA