

Product Data Sheet

Product Name: Activated Protein C (390-404), human TFA

Cat. No.: GC35246

Chemical Properties

Cas. No.

Formula C₉₃H₁₃₁F₃N₂₂O₂₅ M.Wt 2014.16

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Activated Protein C (390-404), human TFA 氨基酸 K 氨基酸 C (APC) 氨基酸 APC 氨基酸[1]

[1]. Mesters RM, et al. Identification of a sequence of human activated protein C (residues 390-404) essential for its anticoagulant activity. J Biol Chem. 1991 Dec 25;266(36):24514-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA